Mediblink[®]
medical devices

EN Mediblink Non-Contact Infrared Body Thermometer M340

INSTRUCTIONS FOR USE

Please read the instructions for use carefully before using the product

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Pažljivo pročitajte upute prije upotrebe uređaja

TABLE OF CONTENTS

Safety Information	6
1. Overview	9
2. Operation	15
3. Troubleshooting	20
4. Replacing the Battery	21
5. Cleaning, Care and Storage	22
6. Disposal	23
7. Warranty	24
8. EMC Declaration	26
9. Body Temperature	30

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Foreword

The Non-Contact Infrared Body Thermometer operating Instructions intend to provide the necessary information for proper operation of M340, thermometer model HTD8816C.

Only body mode was reviewed and certified by notified body.

General knowledge of infrared thermometer and an understanding of the features and functions of the HTD8816C thermometer model are prerequisites for proper use.

The Non-Contact Infrared Body Thermometer is a medical device, and can be used repeatedly, which using life is 5 years.

Please read the manual first before using it, if not fully understand the usages, please stop using the thermometer.

Do not operate any of HTD8816C (REF: M340) thermometer without completely reading and understanding these instructions.

Notice

Purchase or possession of this device does not carry any express or implied license to use with replacement parts which would, alone or in combination with this device, fall within the scope of one of the relating patents.

Safety Information

This device may only be used for the purposes described in these instructions. The manufacturer cannot be held liable for damage caused by incorrect application.

The Non-Contact Infrared Body Thermometer is designed to minimize the possibility of hazards from errors in the software program by following sound and light engineering design processes, Risk Analysis and Software Validation.

Warning

Warnings are identified by the WARNING symbol shown above.

- The Non-Contact Infrared Body Thermometer is to be operated by consumers in the home setting and primary care setting as screening tool. This manual, accessories, Directions for Use, all precautionary information, and specifications should be read before use.
- This product is designed to measure human body temperature on the forehead. Do not use it for any other purpose.
- This product is intended in the home setting and primary care setting as screening tool.
- Do not use the thermometer if it malfunctions or has been damaged in any matter.
- When the ambient temperature of the thermometer changes too much, such as moving the thermometer from one place of lower temperature to another place of higher temperature,

allow the thermometer to remain in a room for 30 minutes where the temperature is between 15°C to 40°C (59–104°F).

- Remove primary batteries if equipment is not likely to be used for long time.
- This product is not waterproof, do not be immersed in water or other liquid; If cleaning and disinfection, please follow the “Care and Storage” section requirements.
- Do not touch the sensor of infrared detection with your fingers.
- If a cold compress on the forehead with patients with fever is used, or other measures to cool down the temperature are taken, measuring data will be lower as expected, and measuring body temperature with this method should be avoided.
- If measuring human forehead temperature, please select “body” mode; for measure other objects, liquids, food and other temperature please select “surface” mode.
- This product must be operated in a stable environment. If the ambient environment is not stable, please check whether there is fog on the sensor. If there is any fog, remove it in accordance with “Care and Storage” section of this instructions for use.
- To avoid the impact on the accuracy of the measurement data, do not use the thermometer near strong electrostatic field or strong magnetic fields.
- Do not mix the old and new batteries to avoid damage to the product.
- It may affect the accuracy of measurements when the forehead is covered by hair, perspiration, cap or scarf.
- The measuring result of this product is only for your reference. If you have any doubt, please measure the temperature using other methods.

 Warning

- The device should be kept out of the reach of children/pets. When not in use, store the device in a dry room and protect it against extreme moisture, heat, lint, dust and direct sunlight. Never place any heavy objects on the storage case.
- Do not throw batteries into fire.
- Use only recommended batteries. Do not use rechargeable batteries.
- This thermometer is not irreplaceable for diagnostic in hospitals.
- Do not fall, disassemble or modify the device.
- Do not use this device if you think it is damaged or if you notice anything unusual.
- This device comprises sensitive components and must be treated with caution. Observe the storage and operating conditions described in the 'Technical Specifications' section.
- Do not perform service/ maintenance while the thermometer is in use.
- When the thermometer is in use, the patient and batteries must not be touched simultaneously.
- Do not use the device if it is damaged/ degraded/ loosened in any way. The continuous use of a damaged unit may cause injury, improper results, or serious danger.
- Based on the current science and technology, other potential allergic reactions are unknown.
- This equipment needs to be installed and put into service in accordance with the information provided in the ACCOMPANYING DOCUMENTS.

1. Overview

Intended Use

The Mediblink Non-Contact Infrared Body Thermometers are designed to be used for intermittent measurement and monitoring of human body temperature by consumers in the home setting and primary care setting as screening tool.

Description of Non-Contact Infrared Body Thermometer

• **Device principle and introduction**

The Mediblink Non-Contact Infrared Body Thermometers are hand-held, reusable, battery operated devices, which can measure human body temperature on forehead, the skin temperature on one's forehead.

The operation principle is based on Infrared Sensor technology. The IR sensor can output different signal when measuring different object temperature or in different ambient temperature, and the ASIC can turn the signal from IR Sensor to a digital value and display it on the LCD.

Description on Controls, Indicators, and Symbols


Figure 1: Forehead thermometer M340

- | | |
|---------------------------------|-------------------------------------|
| 1. Liquid crystal display (LCD) | 10. Indicator of measurement result |
| 2. Battery Cover | 11. Low Battery indicator |
| 3. On/Scan button | 12. Voice on /off indicator |
| 4. MODE button | 13. Fahrenheit |
| 5. MEMO button | 14. Memory number |
| 6. IR sensor | 15. Celsius |
| 7. Surface mode | 16. Memory indicator |
| 8. Body mode | |
| 9. Data indicator | |

Thermometer Applications

Thermometer Model Number	Thermometer Style	Adult		Pediatric	
		Ear	Forehead	Ear	Forehead
REF: M340 Model: HTD8816C	Non-Contact Infrared Body Thermometer		✓		✓

Remark

	Warning		Restriction of Hazardous Substances
	Non-sterile		Batch number*
	Before use, consult instructions for use or electronic instructions for use: http://www.mediblink.com/f/m340.pdf		Compliance with WEEE Standard
	Temperature limit		Atmospheric pressure limitation
	Humidity limitation		CE Marking of Conformity
	This device complies with Part 15 of FCC Rules		Authorized Representative in EU
	Manufacturer		Serial number
	Number of products in one packaging		Medical device
IP22	Degrees of protection: IP22: The first number 2: Protected against solid foreign objects of Φ 12.5 mm and greater. The second number: Protected against vertically falling water drops when enclosure tilted up to 15°.		

*Date of production: the first two digits of the LOT number represent the Year, second two digits Month, and the third two digits represent the Day of production. Example: LOT A221016 = 16th October 2022

Technical Specifications

Measurement Unit	°C / °F
Operating mode	Adjusted mode (Body mode) Direct mode (Surface mode)
Reference Body Site	Axillary
Measuring site	Forehead
Rated output range	Body mode: 34.0–43.0°C / 93.2–109.4°F Surface mode: 0–100.0°C / 32–212°F
Output Range	Body mode: 34.0–43.0°C / 93.2–109.4°F Surface mode: 0–100.0°C / 32–212 °F
Laboratory Accuracy	Body mode: 34.0–34.9°C: ±0.3°C / 93.2–94.8°F: ±0.5°F; 35.0–42.0°C: ±0.2°C / 95.0–107.6°F: ±0.4°F; 42.1–43.0°C: ±0.3°C / 107.8–109.4°F: ±0.5°F Surface mode: ±2°C / ±3.6°F
Display Resolution	0.1°C / 0.1°F
Three-color Backlight (Color Alarm)	35.5–37.3°C / 95.9–99.1°F: Green (Normal Temperature); 37.4–38.0°C (Alarm point) / 99.3–100.4°F: Yellow (Slight Fever) 38.1–43.0°C / 100.6–109.4°F: Red (High Fever) Note: 1. Surface mode is always with Green backlight. 2. In body mode 34.0–35.4°C is with green backlight
Auto Power Off Time	≤18s
Measuring Time	≤2s
Measuring Distance	1–5cm (0,4–2,0in)
Memories	50

Power Supply Requirements	
Batteries	1.5V (AAA) Alkaline battery X 2 (IEC Type LR03)
Adaptable Range	2.6–3.6V
Environmental	
Operating Condition	Operating Temperature: 15°C–40°C (59°F–104°F), Relative Humidity ≤85%, atmospheric pressure: 70Kpa–106Kpa
Transport and Storage Condition	Storage Temperature: –20°C–55°C / –4°F–131°F, Relative Humidity ≤93%, atmospheric pressure: 70Kpa–106Kpa

Dimension and Weighting	
Weight (without batteries)	80g
Size	L:122mm × W:96 mm × H:42mm
Compliance	
Item	Compliant with
Equipment classification	Safety Standards: EN 60601-1: 2006+A1:2013, EN 60601-1-2: 2015
Type of protection	Internally powered equipment (on battery power)
Degree of protection	Non Applied part
Front panel and case labeling	EN ISO15223-1:2016
Temperature	EN ISO80601-2-56:2017
Home healthcare environment	EN 60601-1-11:2015

Calculated values of the indicators according to ISO 80601-2-56

Forehead mode:

	Group A1	Group A2	Group B	Group C
Bias	0.015	0.000	-0.042	-0.040
Standard Deviation	0.123	0.121	0.121	0.137
Limits of Agreement	0.245	0.243	0.243	0.274
Clinical Repeatability (for All Group)	0.071			

Note: the above value is calculated from clinical data of HTD8818A.

Safety classification of ME equipment

Protection against electric shock	Internally powered ME equipment
Applied part	Non Applied part
Protection against harmful ingress of water or particulate matter	IP22
Mode of operation	Continuous operation

Note: Not intended to be sterilized. Not for use in an OXYGEN RICH ENVIRONMENT

2. Operation

2.1 Battery installation

Caution: The Non-Contact Infrared Body Thermometer does not operate with dead batteries and does not input outer power. Install new batteries.

- 1) Pull the battery cover downward, toward the bottom of the Non-Contact Infrared Body Thermometer, and remove the battery through access door;
- 2) Insert two pieces AAA size batteries according to the “+” and “-”;
- 3) Close the battery cover.

2.2 How to Operate

Before Applying the thermometer

Be sure to read and understand all warnings listed of the instructions before use.

- Align the thermometer with the middle of the forehead to measure body temperature (between the eyebrows) and keep the thermometer in vertical position, at a distance of 1–5 cm from the forehead, then press the On/Scan button, and body temperature will be displayed in less then 2 seconds (figure 2).


Figure 2: Measuring position and distance

- When the ambient temperature of the thermometer changes too much, such as when moving the thermometer from one place of lower temperature to another place of higher temperature, allow the thermometer to remain in a room for 30 minutes where the temperature is between 15°C to 40°C.
- The ambient temperature around the test person should be stable. Thermometer should be kept away from the larger flow fan, air-conditioning vents and so on.
- When people are moving from one place of lower temperature to another place of higher temperature, they should remain in the test environment more than 5 minutes, to be consistent with the ambient temperature after the re-measurement.
- Wait at least 1 second for the next measurement. If five continuous measurements are performed it is recommended to wait at least 30 seconds before a new measurement is made.
- You cannot use the thermometer in place where the sun is strong.
- If for some reason the low forehead temperature measurement does not work, try to align the thermometer behind the ears. See figure 3.


Figure 3: Alignment behind the ears

General Setup and Use

• Start measuring

1. Turn on the thermometer by pressing the On/Scan button. The thermometer will perform self-test with all segments displayed for 2 seconds.
2. Align the thermometer to the forehead to keep the distance 1–5 cm, and then press the On/Scan button to start the measurement, read the data.

Note:

- 1) After all symbols on the LCD screen disappear, you will hear a rattle or “beep beep beep” three times, which means that the thermometer self-test has been completed, and the target value of the measured temperature will be displayed on the LCD, while backlight will be displayed according to the appropriate setting among one of the three colors red, green, yellow. Green backlight means ready for next measurement. When 37.4–38.0°C, backlight is yellow, meaning slight fever warning. Please pay attention to body temperature. When the body temperature is above 38.1°C, backlight is red, meaning high fever. Please take action to cool down or go to see a doctor.

2) To ensure the accuracy of the measurement, wait at least 30 seconds after 5 consecutive measurements.

- **Mode conversion**

When the device is ON press the MODE button to make conversion between “body” mode and “surface” mode.

“Body” mode is used to measure human body temperature, the “surface” mode is used to measure the surface temperature (the factory default is “body” mode).

- **Recalling and Erasing Memory Data**

The last temperature taken before the thermometer powers off is stored in memory, and up to 50 measurements are stored in memory.

Steps:

- 1) In the boot or shutdown state, short press the MEMO button to view the history of measured values.
- 2) An empty memory cell shows “---°C” or “---°F”.
- 3) Temperature readings can be stored in memory. Up to 50 temperature readings can be stored into the memory cells and new readings automatically overwrite historical data.
- 4) In boot mode, long-press the MEMO button until the LCD displays “CLR”, which means that all stored data is erased completely.

- **Parameter settings**

Steps:

- 1) Unit Set – F1

Under the boot mode, long press the MODE button to enter F1, press the “MODE” or “MEMO” button to switch

between Celsius and Fahrenheit temperature units, press the On/Scan button to confirm the unit settings (factory default is Celsius).

2) Fever alert Set – F2

Under F1 state, press On/Scan button to enter the F2, press the “MODE” button to plus 0.1°C, press the “MEMO” button decrease 0.1°C, long press to accelerate the speed of temperature regulation, and finally press the On/Scan button to save the settings (the factory default is 38.1°C).

3) Prompt sound settings – F3

Under F2 state, press On/Scan button to enter F3, press MODE button or MEMO button to disable or enable the voice, and press the On/Scan button to confirm the settings (the factory default is enabled voice).


• **Prompt sound settings ON/OFF function**


Under the boot mode, short press power on/off button, and at the same time press the MODE button to take on or take off the Prompt sound function.

• **Restore to factory default settings function**

Under the boot mode, short press power on/off button, at the same time press the MEMO button. Until LCD display “rst”. Two seconds later, F1-F3 parameters will be back at factory default settings.

3. Troubleshooting

Message	Situation	Solution
Body °C 	Temperature taken is not within typical human temperature range (34.0~43.0°C or 93.2~109.4°F).	Make sure that temperature is measured in the forehead, and not on other human body site.
Body °C 	Measured over the distance: 1–5 cm (0.4–2 in).	Optimum measurement distance is 1 cm.
Body °C 	Incorrect test position.	See figure 2: Measuring position and distance.
	Subjects such as: forehead hair, antipyretic stickers, head with sweat, etc.	Subjects sit quietly 5–10 minutes before the test.
	Some people's body temperature is lower than the general population.	The main concern is fever – high temperature.
Body °C 	Operating temperature exceeds the range of specified temperature.	Move to a room within the operating range and wait 30 minutes before taking temperature.
	The screen flickers and automatically turns off.	Replace the battery or the product has been damaged, and needs a repair.

	Battery capacity is too low. Taking temperature is not allowed.	Insert a new battery.
	Ambient temperature changes too fast.	Wait until the ambient temperature is stable.
	(1) Power is off. (2) Improper battery installation. (3) The battery is exhausted. (4) Display remains blank.	(1) Press ON button again. (2) Check the battery polarity. (3) Replace with a new battery. (4) Contact the retailer or service center.

4. Replacing the Battery

1. Open and release battery cover following indicator on the surface of battery cover. Before changing the battery be sure the system is already powered off.
2. Remove the battery and replace them with 2 new ones, type AAA, Make sure to align the batteries properly as it is indicated inside the battery cover.
3. Slide the battery cover back in, until it snaps in place.
Do not dispose used batteries in household waste. Take them to special local collection sites.
4. In case, if system is latched up after changing battery. You may not follow up the process of rule one. Just take off battery, wait for 30 sec, then load battery again.

**Do not recharge, disassemble or dispose of in fire.**

1. The typical service life of the new and unused batteries is 2000 measurements if the operation time is 18s.
2. Use only recommended batteries, do not recharge non-rechargeable batteries and do not burn them.
3. Remove the batteries if the thermometer is not going to be used for a longer period.

5. Cleaning, Care and Storage

The lens are very delicate.

It is very important to protect the lens from dirt and damage.

Use a clean, soft cloth to clean the surface of the device and LCD. Do not use solvents nor immerse the device into water or other liquids.

Always keep the thermometer within the storage temperature range (-20°C to 55°C or -4°F to 131°F) and humidity range ($\leq 93\%$ non-condensing).

It is recommended to store the thermometer in a dry location free from dust. Do not expose the thermometer to direct sunlight, high temperature/ humidity or any extreme environment, otherwise the function will be reduced.

When the ambient temperature of the thermometer changes too much, such as moving the thermometer from one place of

lower temperature to another place of higher temperature, allow the thermometer to remain in a room for 30 minutes where the temperature is between 15°C to 40°C.

6. Disposal

- Used batteries should not be disposed of in the household rubbish. Used batteries should be deposited at a collection point.
- At the end of its life, the appliance should not be disposed in household rubbish. Enquire about the options for environment-friendly and appropriate disposal. Take local regulations into account.

7. Warranty

Product: Mediblink Non-Contact Infrared Body Thermometer M340

Manufactured for (EU importer and distributor): Mediblink d.o.o., Gubčeva cesta 19, 8210 Trebnje, Slovenia; info@mediblink.com; www.mediblink.com

Sellers name, address, signature and stamp*:

Date of extradition/sales*:

*If the invoice is accompanied by this warranty, and if all above information can be seen from the invoice, it is not necessary to fill in this field.

WARRANTY TERMS

Dear customers!

The warranty period is **5 years** and starts from the day of product purchase. In case of product claim, you have to show the invoice. We kindly ask you to save the invoice!

Unfortunately, wrong handling with the device is a reason for 95% of customer complains. You can easily avoid any problem, by getting useful information provided by our special service department. To reach our service department, you can call or send e-mail to Mediblink local distributor.

Before sending the product back to retailer, we kindly ask you to call our service department, to get help about how to use the device to save you with unneeded trips.

The manufacturer guarantees free elimination of all imperfections due to defects in material or manufacturing procedure by repairing or replacing the product. In case that the product can not be repaired or replaced, the customer will get the money refund. The guarantee is not valid in case of the force majeure, accidents or unexpected events (such as lightning, water, fire etc.), incorrect use or incorrect transport, non-compliance with safety and maintaining regulations or in case of unprofessional product intervention.

Traces of every day product usage (scratches, abrasions) and not subject to claim. The warranty does not eliminate the customer rights, which originate from seller responsibility for product flaws. By accepting the claimed product by the service

department, the service department does not take responsibility for loss of saved data or settings on the product. All product repairs, which are performed out of product warranty period, have to be paid by customer by prior notice.

The manufacturer guarantees the product quality and flawless product operation in the warranty period, which starts with the day of product purchase. If the product can not be repaired in 45 days, the product will be replaced with a new one. In case that the product can not be replaced, the money will be refunded to the customer.

In case of product claim, call or send e-mail to Mediblink local distributor.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

The thermometer is initially calibrated at the time of manufacture. If this thermometer is used according to the instructions for use, a periodic re-adjustment is not required. If at any time you question the accuracy of temperature measurements, please contact the manufacturer.

Manufacturer:
HeTaiDa Technology Co., Ltd
Room 801, 802, 803, 804, 901, 2#
Building Scientific Research Center,
Songhu Intelligent Valley, No.6 Minfu
Road, Liaobu Town, Dongguan City,
Guangdong Province, P.R.China

Manufactured for
(EU importer and distributor):
Mediblink d.o.o.
Gubčeva cesta 19
8210 Trebnje
Slovenia
info@mediblink.com
www.mediblink.com

Share Info Consultant Service
LLC Repräsentanzbüro
Heerdtter Lohweg 83,
40549 Düsseldorf,
Germany

Software version No.: V31
Instructions for use version number: A.10
Issue date: 07. 04. 2021
Date of last modification: 16. 10. 2023

3.0, 16. 10. 2023

8. EMC Declaration

- 1) This equipment needs to be installed and put into service in accordance with the information provided in the ACCOMPANYING DOCUMENTS.

This product needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided, and this unit can be affected by portable and mobile RF communications equipment.

- 2)* Caution: Do not use a mobile phone or other devices that emit electromagnetic fields, near the unit. This may result in incorrect operation of the unit.
- 3)* Caution: This unit has been thoroughly tested and inspected to assure proper performance and operation!
- 4)* Caution: This machine should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this machine should be observed to verify normal operation in the configuration in which it will be used.

Guidance and manufacture's declaration – electromagnetic emission

The infrared body thermometer is intended for use in the electromagnetic environment specified below. The customer of the user of the infrared body thermometer should assure that it is used in such an environment.

Emission test	Compliance
RF emissions CISPR 11	Group 1
RF emission CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	Not applicable
Voltage fluctuations / flicker emissions IEC 61000-3-3	Not applicable

Guidance and manufacture's declaration – electromagnetic immunity

The infrared body thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of infrared body thermometer should assure that it is used in such an environment.

Anti-interference detection	IEC 60601 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	Contact: ± 8 KV Air: $\pm 2, \pm 4, \pm 8, \pm 15$ KV	Contact: ± 8 KV Air: $\pm 2, \pm 4, \pm 8, \pm 15$ KV
Electrical fast transient / burst IEC 61000-4-4	The input a.c. power ports: ± 2 KV The input d.c. power ports: ± 2 KV Signal input/output ports: ± 1 KV	Not applicable
Surge IEC 61000-4-5	Input power ports: +0.5, +1.0 KV Signal input/output: +2.0 KV	Not applicable
Voltage dips IEC 61000-4-11	0.5 cycles for >95% (sync angle (degrees): 0, 45, 90, 135, 180, 225, 270, 315) 1 cycles for >95% UT (sync angle (degrees): 0) 25 (50Hz) / 30 (60Hz) cycles for 30% U T (sync angle (degrees): 0)	Not applicable
Voltage interruption IEC 61000-4-11	250 (50Hz) / 300 (60Hz) cycles for >95% UT (sync angle (degrees):0)	
Power frequency (50Hz / 60Hz) magnetic field IEC 61000-4-8	30A/m	30A/m

NOTE: UT is the a.c. mains voltage prior to application of the test level.

Guidance and manufacture's declaration – electromagnetic immunity

The infrared body thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of the infrared body thermometer should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level
Conducted RF IEC 61000-4-6	3Vrms 150kHz to 80MHz	Not applicable
Radiated RF IEC 61000-4-3	Professional healthcare environment: 3V/m Home healthcare environment: 10Vm 80MHz to 2700MHz	Professional healthcare environment: 3V/m Home healthcare environment: 10 Vm 80 MHz to 2700MHz (G)

NOTE 1: At 80MHz and 800MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Non-Contact Infrared Body Thermometer is used exceeds the applicable RF compliance level above, the Non-Contact Infrared Body Thermometer should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Non-Contact Infrared Body Thermometer.
- b) Over the frequency range 150kHz to 80MHz, field strengths should be less than 3V/m.

Guidance and manufacturer's declaration – RF wireless communication equipment immunity

Test frequency (MHz)	Band ^a (MHz)	Service ^a	Modulation ^b	Maximum power (W)	Distance (m)	Immunity test level (V/m)
385	380–390	TETRA 400	Pulse modulation ^b 18 Hz	1,8	0,3	27
450	430–470	GMRS 460, FRS 460	FM ^c ± 5 kHz deviation 1 kHz sine	2	0,3	28
710 745 780	704–787	LTE Band 13, 17	Pulse modulation ^b 217 Hz	0,2	0,3	9
810 870 930	800–960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation ^b 18 Hz	2	0,3	28
1720 1845 1970	1700–1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation ^b 217 Hz	2	0,3	28
2450	2400– 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ^b 217 Hz	2	0,3	28
5240 5500 5785	5100– 5800	WLAN 802.11 a/n	Pulse modulation ^b 217 Hz	0,2	0,3	9

^a For some services, only the uplink frequencies are included.

^b The carrier shall be modulated using a 50 % duty cycle square wave signal.

^c As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

9. Body Temperature

- Body temperature varies from person to person and fluctuates during the course of the day. For this reason, it is suggested to know one's normal, healthy forehead temperature to correctly determine the temperature.
- Body temperature runs approximately from 35.5°C to 37.8°C (95.9°F–100°F). To determine if someone has a fever, compare the temperature detected with the person's normal temperature. A rise over the reference body temperature of 1°C (1°F) or more is generally indication of fever.
- Different measurement sites (rectal, axillary, oral, frontal, auricular) will give different readings. Therefore it is wrong to compare the measurement taken from different sites.
- Here below are typical temperatures for adults, based on different measurement sites:
 - Rectal: 36.6 to 38°C / 97.9–100.4°F
 - Axillary: 34.7 to 37.3°C / 94.5–99.1°F

ASTM laboratory accuracy requirements in the display range of 37 to 38.9°C (98 to 102°F) for IR thermometers is $\pm 0.2^{\circ}\text{C}$ ($\pm 0.4^{\circ}\text{F}$), whereas for mercury in-glass and electronic thermometers, the requirement per ASTM Standards E667-86 and E1112-86 is $\pm 0.1^{\circ}\text{C}$ ($\pm 0.2^{\circ}\text{F}$).

Caution: This infrared thermometer meets requirements established in ASTM Standard (E1965-98), except of clause 5.2.2. It displays subject's temperature over a range of 34.0~43.0°C. Full responsibility for the conformance of this product to the standard is assumed by Hetaida Technology Co., (Ltd. Add:

Room 801, 802, 803, 804, 901, 2# Building Scientific Research Center, Songhu Intelligent Valley, No.6 Minfu Road, Liaobu Town, Dongguan City, Guangdong Province, P.R.China).

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