Mediblink®

medical devices

EN

Mediblink Digital thermometer with flexible tip Flexy M374

INSTRUCTIONS FOR USE

PLEASE READ CAREFULLY BEFORE USE

Mediblink Digital thermometer with flexible tip Flexy M374 provides a quick and highly accurate reading of an individual's body temperature. To better understand its functions and to provide years of dependable results, please read all instructions for use first.

Content: 1 Digital thermometer, 1 User's Manual, 1 Storage Case

INTENDED USE

The digital thermometers are intended to measure the human body temperature in regular mode orally, rectally, under the arm. The devices are reusable for clinical or home use for people of all ages, including children under 8 years old with adult supervision.

INTENDED USER/INDICATION FOR USE

Patients or healthy people who want to measure their body temperature, and medical workers who want to measure the patient's body temperature.

Patient target group: Patients or healthy people who need or want to have their temperature measured.

PRODUCT ILLUSTRATION


Figure 1: Digital thermometer components

SPECIFICATIONS

Product name:

Mediblink Digital thermometer with flexible tip Flexy M374

Type:

Digital thermometer (Not Predictive)

Reference number / Model:

M374 / DMT-4362

Measure Range:

32.0°C–42.9°C (90.0°F–109.9°F) (°C/°F chosen by manufacturer)

Accuracy:

±0.1°C (±0.2°F) during 35.5°C–42.0°C (95.9°F–107.6°F) at 18°C–28°C (64.4°F–82.4°F) ambient operating range, ±0.2°C (±0.4°F) for other measuring and ambient operating range

Operating mode:

Direct Mode

Display:

Liquid crystal display, 3 1/2 digits

Memory:

For storing the last measured value

Battery:

One 1.5V DC button battery (size LR41 or SR41, UCC 392)

Battery life:

Approx. 200 hours of continuous operation or 1 year with 3 measurements per day

Dimension:

14.3 x 2.5 x 1.4 cm (L x W x H)

Weight:

Approx. 16 grams including battery

Expected service life:

Three years

Ambient operating range:

Temperature: 5°C–40°C (41°F–104°F) Relative humidity: 15%–95% RH Atmospheric Pressure: 700 hPa–1060 hPa

Storage and transportation condition:

Temperature: -20°C–55°C (-4°F–131°F) Relative humidity: 15%–95% RH Atmospheric Pressure: 700 hPa–1060 hPa

Ingress Protection Rating:

IP 27

Classification:

Type BF

Contraindications:

No contraindications

Conforming to standards:

ISO 80601-2-56 Medical electrical equipment – Part 2-56: Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement, EN 60601-1-11 Medical electrical equipment – Part 1-11: General requirements for basic safety and essential performance – Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment and complies with the requirements of EN 60601-1-2 (EMC), IEC/EN 60601-1 (Safety) standards. And the manufacturer is ISO 13485 certified.

⚠ PRECAUTIONS

1. The performance of the device may be degraded should one or more of the following occur:

- Operation outside the manufacturer's stated temperature and humidity range.
- Storage outside the manufacturer's stated temperature and humidity range.
- Mechanical shock (for example, drop test) or degraded sensor.
- Patient temperature is below ambient temperature.

2. Portable and mobile RF communications can affect the device. The device needs special pre-cautions regarding EMC according to the EMC information provided in the accompany documents.

⚠ WARNINGS

- Read instructions thoroughly before using digital thermometer.
- Choking Hazard: Thermometer cover and battery may be fatal if swallowed. Do not allow children to use this device without parental supervision.
- Do not use thermometer in ear. Designed use is for oral, rectal, and armpit (axilla) readings only.
- Do not place thermometer battery near extreme heat as it may explode.
- Remove battery from the device when not in operation for a long time.
- The use of temperature readings for self-diagnosis is dangerous. Consult your doctor for the interpretation of results. Self-diagnosis may lead to the worsening of existing disease conditions.
- Do not attempt measurements when the thermometer is wet as inaccurate readings may result.
- Do not bite the thermometer. Doing so may lead to leakage and/or injury.
- Do not attempt to disassemble or repair the thermometer. Doing so may result in inaccurate readings.
- After each use, disinfect the thermometer especially in case the device is used by more than one person.
- Do not force the thermometer into the rectum. Stop insertion and abort the measurement when pain is present. Failure to do so may lead to injury.
- Do not use thermometer orally after being used rectally.
- For children who are two years old or younger, please do not use the devices orally.
- If the unit has been stored at temperatures over 5°C–40°C (41°F–104°F), leave it in 5°C–40°C (41°F–104°F) ambient temperature for about 15 minutes before using it.
- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- PORTABLE RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the [ME EQUIPMENT or ME SYSTEM], including cables specified by the MANUFACTURER. Otherwise, degradation of the performance of this equipment could result.
- It is not intended for use in the oxygen rich environment and presence of flammable anaesthetic mixture with air, oxygen or nitrous oxide.
- Do not put the thermometer in direct sunlight or in cotton wool, otherwise the accuracy will be affected.
- ME equipment should not be cleaned and disinfected while in use.

°C/°F SWITCHABLE

Temperature readings are available in the Celsius scale (°C/°F), located in the upper right corner of LCD.) With the unit off, press and hold the On/Off Button for approximately 2 seconds to change the current setting.

DIRECTIONS FOR USE

1. Press the On/Off Button next to LCD display in Figure 1. A tone will sound as the screen shows 88.8°, followed by last recorded temperature. After showing the self-test temperature, the thermometer is now in the testing mode.

2. Position thermometer in desired location (mouth, rectum, or armpit).

a) **Oral Use:** Place thermometer under tongue as indicated by "✓" position shown in Figure 2. Close your mouth and breathe evenly through the nose to prevent the measurement from being influenced by inhaled/exhaled air.


Figure 2

b) **Rectal Use:** Lubricate silver probe tip with petroleum jelly for easy insertion. Gently insert sensor approximately 1cm (less than 1/2") into rectum.

c) **Armpit Use:** Wipe armpit dry. Place probe in armpit and keep arm pressed firmly at side. From a medical viewpoint, this method will always provide inaccurate readings, and should not be used if precise measurements are required.

3. The degree sign flashes throughout the testing process. When flashing stops an alarm will beep for approximately 10 seconds. The measured reading will appear on the LCD simultaneously. The minimum measurement time until the signalling tone (beep) must be maintained without exception. The measurement continues even after the buzzer notification.

IMPORTANT: in order to achieve better body temperature measurement result, it is recommended to keep the probe in mouth and rectum about 2 minutes, or in armpit about 5 minutes regardless of the beep sound and at least 30 seconds measurement interval should be maintained.

*Note: Normally the buzzes are "Bi-Bi-Bi-Bi-"; alarm beeps more rapidly when temperature reaches 37.8°C (100°F) or higher and the buzzers are " Bi-Bi-Bi- Bi- Bi-Bi- Bi-Bi-Bi- Bi-Bi-Bi"

4. To prolong battery life, press the On/Off Button to turn unit off after testing is complete. If no action is taken, the unit will automatically shut off after around 10 minutes.

TROUBLESHOOTING

Error message	Problem	Solution
Lo	Temperature taken is lower than 32.0°C (90.0°F)	Turn off, wait one minute and take a new temperature via close contact and sufficient rest.
Hi	Temperature taken is higher than 42.9°C (109.9°F)	Turn off, wait one minute and take a new temperature via close contact and sufficient rest.
Err	The system is not functioning properly.	Unload the battery, wait for 1 minute and repower it. If the message reappears, contact the retailer for service.
	Dead battery: Battery icon is flashing, can't measure.	Replace the battery.

BATTERY REPLACEMENT

1. Replace battery when the symbol in the lower right corner of the LCD display starts flashing.

2. Pull battery cover off as shown in Figure 3.

3. Gently pull-out plastic circuit board with battery chamber approximately 1cm (slightly less than 1/2") (Figure 4).

4. Use pointed object such as a pen to remove old battery. Discard battery lawfully. Replace with new 1.5V DC button type LR41 or SR41, UCC392, or equivalent. Be sure battery is installed with polarity facing up. (Figure 5).

5. Slide battery chamber back into place and attach cover.

Note: The following schematic diagram of battery replacement is from the back of the thermometer.

Figure 3

Figure 4

Figure 5

CALIBRATION

The thermometer is initially calibrated at the time of manufacture. If the thermometer is used according to the instruction of use, periodic readjustment is not required. However, we recommend checking calibration every two years or whenever clinical accuracy of the thermometer is in question. Set the temperature of the high-precision water bath temperature bath to 33.00°C, 37.00°C, 42.00°C. Turn on the thermometer and insert into the high-precision water bath and then check the laboratory accuracy of thermometer. Repeat the previous measurement process after it is stable for 10 minutes. Please send the complete device to the dealers or manufacturer.

The above recommendations do not supersede the legal requirements. The user must always comply with legal requirements for the control of the measurement, functionality and accuracy of the device which are required by the scope of relevant laws, directives or ordinances where the device is used.

CLEANING AND DISINFECTION

1. Pre-clean the surface of the thermometer probe using potable water or wipe it with a clean, soft cloth dampened with potable water to remove surface dirt.

2. Immerse the thermometer probe completely in high-purity water for at least 1 minute.

3. Gently wipe the thermometer with a dry, clean, soft cloth to remove any residue.

4. Repeat steps 2 and 3 three times until no dirt is visible upon visual inspection after cleaning.

5. Disinfect the probe:

- After use in the armpit, mouth, or rectum:

Method A (High-level disinfection): Immerse the thermometer probe in 0.55% OPA (Ortho-phthalaldehyde), such as CIDEX OPA, for at least 12 minutes at 20°C (68°F).
Use this method before oral or rectal use.

Method B (Intermediate disinfection): Using a clean soft cloth dipped in 70% isopropanol, wipe the probe three times, for at least one minute each time.

6. Repeat steps 2 to 4 to remove any disinfectant residue.

7. After cleaning and disinfection are complete, check that the product appears normal. Then, place the thermometer in its dust cover (if available) and store it in a cool, dry, clean environment, away from high temperatures and corrosive substances.

Precautions:

1. The device must be cleaned and disinfected before and after each use.

2. For oral or rectal measurements, advanced disinfection is strongly recommended.

3. For the disinfectants or cleaning agents mentioned above, always follow the manufacturers' instructions.

4. The water temperature used in any of the above processes should not exceed 45°C (113°F).

5. Do not use benzene, paint thinner, gasoline, or other strong solvents to clean the thermometer.

6. Do not attempt to disinfect the sensing probe (tip) of the thermometer by immersing it in 70% isopropanol, OPA, or hot water (above 45°C / 113°F) for extended periods.

7. Do not use ultrasound or any automatic/semi-automatic cleaning methods on the thermometer.

8. Do not use any disinfection methods other than the recommended ones. Unapproved methods may cause permanent damage to the equipment.

9. During its service life, the product can be reused approximately 3,000 times if the above cleaning and disinfection procedures are followed.

ELECTROMAGNETIC COMPATIBILITY INFORMATION

The device satisfies the EMC requirements of the international standard IEC 60601-1-2. The requirements are satisfied under conditions described in the table below. The device is an electrical medical product and is subject to special precautionary measures with regard to EMC which must be published in the instructions for use. Portable and mobile HF communications equipment can affect the device. Use of the unit in conjunction with non-approved accessories can affect the device negatively and alter the electromagnetic compatibility. The device should not be used directly adjacent to or between other electrical equipment.

Electromagnetic Compatibility Information – electromagnetic emission		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.		
Emission test	Compliance	
Conducted emissions CISPR 11	Not applicable	
RF emissions CISPR 11	Group 1, Class B	
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Not applicable	

Electromagnetic Compatibility Information – electromagnetic immunity		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.		
Immunity test	IEC 60601 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	±8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Electrical transient/burst IEC 61000-4-4	±2 kV for power supply lines 100 kHz repetition frequency ±1 kV for input/output lines	Not applicable
Surge IEC 61000-4-5	+0.5 kV, ±2 kV differential mode line-line	Not applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% U _i (100% dip in U _i) for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% U _i (100% dip in U _i) for 1 cycle at 0° 70% U _i (30% dip in U _i) for 25/30 cycles at 0° 0% U _i (100% dip in U _i) for 250/300 cycles at 0°	Not applicable
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m, 50/60 Hz	30 A/m, 50/60 Hz
NOTE: U _i is the a.c. mains voltage prior to application of the test level.		

Electromagnetic Compatibility Information – electromagnetic immunity		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.		
Immunity test	IEC 60601 test level	Compliance level
Conducted RF IEC 61000-4-6	3 Vrms, 150 kHz to 80 MHz 6 Vrms, 150 kHz to 80 MHz ISM bands	Not applicable
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz	10 V/m
Proximity magnetic fields IEC 61000-4-39	30 kHz, 8 A/m, CW 134.2 kHz, 65 A/m, Pulse modulation 2.1 kHz 13.56 kHz, 7.5 A/m, CW 134.2 kHz, 65 A/m, Pulse modulation 50 kHz	30 kHz, 8 A/m, CW 134.2 kHz, 65 A/m, Pulse modulation 2.1 kHz 13.56 kHz, 7.5 A/m, CW 134.2 kHz, 65 A/m, Pulse modulation 50 kHz

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

⁴ The ISM (industrial, scientific and medical) bands between 0.15 MHz and 80 MHz are 6.765 MHz to 6.795 MHz, 13.553 MHz to 13.567 MHz, 26.957 MHz to 27.283 MHz and 40.66 MHz to 40.70 MHz. The amateur radio bands between 0.15 MHz and 80 MHz are 1.8 MHz to 2.0 MHz, 3.5 MHz to 4.0 MHz, 5.3 MHz to 5.4 MHz, 7.0 MHz to 7.3 MHz, 10.1 MHz to 10.15 MHz, 14 MHz to 14.2 MHz, 18.07 MHz to 18.17 MHz, 21 MHz to 21.4 MHz, 24.89 MHz to 24.99 MHz, 28 MHz to 29.7 MHz and 50 MHz to 54 MHz.

⁵ The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and the frequency range 80 MHz to 2.7 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in these frequency ranges.

⁶ Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radios, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

⁷ Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between RF wireless communications equipment				
The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between RF wireless communications equipment and the device as recommended below, according to the maximum output power of the communications equipment.				
Frequency (MHz)	Maximum Power (W)	Distance (m)	IEC 60601 Test Level	Compliance Level
385	1.8	0.3	27	27
450	2	0.3	28	28
710, 745, 780	0.2	0.3	9	9
810, 870, 930	2	0.3	28	28
1720, 1845, 1970	2	0.3	28	28
2450	2	0.3	28	28
5240, 5500, 5785	0.2	0.3	9	9

Note 1: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- ⚠ WARNING**
- This device should not be used in the vicinity or on the top of other electronic equipment such as cell phone, transceiver or radio control products. If you have to do so, the device should be observed to verify normal operation.
 - The use of accessories and power cord other than those specified, with the exception of cables sold by the manufacturer of the equipment or system as replacement parts for internal components, may result in increased emissions or decreased of the equipment or system.

LIST OF SYMBOLS

Symbol	Explanation
	Electrical devices are recyclable material and should not be disposed of with household waste after use! Help us protect the environment and save resources and take this device to the appropriate collection points. Please contact the organization which is responsible for waste disposal in your area if any questions.
	Applied part of type BF
	MR Unsafe
	Before use, read the instructions. Electronic instructions for use: http://www.mediblink.com/f/m374.pdf
IP27	The first number: 2 - Protected against solid foreign objects of 12.5 mm Ø and greater. The second number: 7 - Protected against the effects of temporary immersion in water.
	Manufacturer
	Authorized Representative in EU
	Warning!
	Batch number*
	Product reference number
	Model number
	Medical device
	Unique Device Identifier
	Atmospheric pressure limitation for Storage and Transportation
	Humidity limitation for Storage and Transportation
	Temperature limitation for Storage and Transportation
	Number of products included in packaging
	CE Marking of Conformity, Notified body ID number

*The production date can be read from LOT number [YYYYMM]; first four digits represent the Year and last two digits Month of production. Example: LOT 202503 = March 2025

WARRANTY

Product: Mediblink Digital thermometer with flexible tip Flexy M374
Manufactured for (importer for EU & distributor): Mediblink d.o.o., Gubčeva cesta 19, 8210 Trebnje, Slovenia; info@mediblink.com; www.mediblink.com

Seller's name, address, signature and stamp*:

Date of extradition / sales*:

*If the invoice is accompanied by this warranty, and if all above information can be seen from the invoice, it is not necessary to fill in this field.

WARRANTY TERMS

Dear customers!

The warranty period is **10 years** and starts from the day of product purchase. In case of product claim, you have to show the invoice. We kindly ask you to save the invoice!

Unfortunately, wrong handling with the device is a reason for 95% of customer complains. You can easily avoid any problem, by getting useful information provided by our special service department. To reach our service department, you can call or send e-mail to Mediblink local distributor.

Before sending the product back to retailer, we kindly ask you to call our service department, to get help about how to use the device to save you with unneeded trips.

The manufacturer guarantees free elimination of all imperfections due to defects in material or manufacturing procedure by repairing or replacing the product. In case that the product cannot be repaired or replaced, the customer will get the money refund. The guarantee is not valid in case of the force majeure, accidents or unexpected events (such as lightning, water, fire etc.), incorrect use or incorrect transport, non-compliance with safety and maintaining regulations or in case of unprofessional product intervention.

Traces of every day product usage (scratches, abrasions) and not subject to claim. The warranty does not eliminate the customer rights, which originate from seller responsibility for product flaws. By accepting the claimed product by the service department, the service department does not take responsibility for loss of saved data or settings on the product. All product repairs, which are performed out of product warranty period, have to be paid by customer by prior notice.

The manufacturer guarantees the product quality and flawless product operation in the warranty period, which starts with the day of product purchase. If the product cannot be repaired in 45 days, the product will be replaced with a new one. In case that the product cannot be replaced, the money will be refunded to the customer.

In case of product claim, call or send e-mail to Mediblink local distributor.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

JOYTECH Healthcare Co., Ltd.
No.365, Wuzhou Road, 311100 Hangzhou, Zhejiang Province, PEOPLE'S REPUBLIC OF CHINA
E-mail: info@sejoy.com

Shanghai International Holding Corp. GmbH (Europe)
Eiffestrasse 80, 20537 Hamburg, Germany
E-mail: shholding@hotmail.com

Manufactured for (Importer for EU & Distributor):
Mediblink d.o.o., Gubčeva cesta 19, 8210 Trebnje, Slovenia,
info@mediblink.com, www.mediblink.com

Instructions for use version: V.01
Issue date: 15. 7. 2021
Date of last modification: 30. 5. 2025

0123

The product is in compliance with the requirements of MDR(EU)2017/745.

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JOYTECH Healthcare Co., Ltd., No.365, Wuzhou Road, 311100 Hangzhou, Zhejiang Province, PEOPLE'S REPUBLIC OF CHINA