Mediblink®

medical devices

EN MEDIBLINK Compressor nebulizer COMPACT M440

INSTRUCTIONS FOR USE

These instructions are important
Please keep them for future reference

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Introduction	6
Device Components	7
Accessories	8
General Warnings	10
2-stage nebuliser cup	11
2 Masks: for adults and children	13
Mouthpiece	13
Preparing and using the appliance	13
Cleaning and maintaining	14
Checking and replacing filters	15
Troubleshooting guide	
Product Specification	16
EMC Information	17
Legend of symbols	21
Conditions for storage and transportation	21
Environmental conditions during use	21
Warranty	22

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Introduction

Dear Customer,

Thank you for choosing Mediblink Compact M440 nebulizer, a piston aerosol therapy unit designed and manufactured according the most up-to-date technologies.

Mediblink Compact M440 is a class IIa medical device which converts medicines from liquid or suspension form to aerosol form, so they can be administered directly to the airways.

Mediblink Compact M440 nebulizer allows for deeper penetration of the medicine in the respiratory tract (MMAD: 2,74 μm), ensured by the perfect size of the aerosol particles.

Warning!


This product must be assembled, maintained and handled by an adult. Never leave children alone during use: the unit contains small parts that could be swallowed.

This product must remain in place for use only for the time strictly needed. Once treatment has been completed, the unit must be stored in a safe place out of reach of children. Before using the appliance, we highly recommend that you read the short list of warnings contained in this manual, to make sure you have correctly understood how to use it.

Device Components

Picture 1: Nebulizer Components


- A – On/Off switch
- B – Compressed air outlet
- C – Filter
- D – Nebuliser cup support
- E – Compressor (main unit)
- F – Openings for cooling air
- G – Power cord

Accessories

Picture 2: Nebulizer Accessories

- H – Adult mask
- I – Children & Baby mask
- J – Mouthpiece
- K – 2-stage Nebulizer cup
- L – Air tube
- M – Filters

General Warnings

1. Before use, carefully read the information in this manual, and keep for future reference.
2. Unit designed for aerosol therapy. Follow the instructions of a medical practitioner as to the type of medicine to use, the dose, frequency and duration of inhalations. Any use other than the one for which the unit is intended is considered improper and therefore dangerous; Mediblink cannot be held liable for any damage caused by incorrect or unreasonable use or if the device is used in electrical systems that do not conform to safety regulations in effect.
3. A power outage, a sudden failure or any other adverse condition may cause the unit to stop working, therefore it is recommended that users have an alternative unit or medicine (agreed upon with a physician) that can be used in the event this occurs.
4. After removing the appliance from its packaging, check that it is intact, without any visible damage that may have occurred during transport. If in doubt, do not use the appliance, and contact the Mediblink representative company in your country.
5. Packaging (bags, box, etc.) may be dangerous and should be kept out of reach of children.
6. Before connecting the appliance, check that your mains voltage corresponds to that shown on the rating label. The rating label is located on the bottom of the appliance.
7. In case of incompatibility between the appliance plug and the electrical socket, use certified adapters if the existing legislation in the individual country allows for it or have an Mediblink authorised technician replace the plug.
8. To prevent overheating and damage to the compressor, switch off the appliance for at least 30 minutes after every 30 minutes of uninterrupted use.
9. Caution! Using the appliance before 30 minutes have elapsed may result in motor overheating and cause the safety cut-off devices to activate.
10. To ensure correct operation, make sure that the air filter is dry.
11. Do not use the appliance in the presence of nitrous oxide, oxygen or flammable air/anaesthetic mixtures.
12. Turn off the appliance and disconnect from the mains after use and before adding more medicine. Do not exceed the maximum level shown on the nebuliser cup when filling.
13. Keep the appliance and power cord away from hot surfaces.
14. Do not use the appliance while having a bath or shower, or when in a damp place, or close to bath-tubs, sinks, washbasins, or in any other situation in which there are liquids that may come into contact with the appliance.
15. Never touch the appliance with wet or damp hands.

General Warnings

16. Do not drop or put the device into water or other liquids. Should this occur, unplug it immediately, stop using the appliance and consult a qualified technician.
17. Do not block the openings for cooling air.
18. Do not use the device if you feel drowsy.
19. If children or persons with physical or mental disabilities use the appliance, proper supervision is required. Children must be kept under adult supervision to ensure that they do not play with the appliance. The appliance contains small parts that can be swallowed and the power cord poses risk of strangling.
20. Do not use connectors or accessories not recommended by the manufacturer.
21. Make sure the appliance is placed on a level, stable surface when used in order to prevent spillage.
22. During use, the surface must be free of any objects that may obstruct the correct flow of air.
23. Use the appliance only in vertical position.
24. Never leave the appliance plugged in when it is not in use or when it is unattended.
25. Before cleaning or performing any maintenance work on the appliance, disconnect it from the mains by pulling the plug out of the electrical socket.
26. If you decide to no longer use the appliance, disconnect from the mains and cut off the power cord so the appliance cannot work. Dispose of the appliance and power cord in compliance with existing laws. We also recommend eliminating any parts that are potential hazards, particularly for children.
27. Refer to laws in force for information regarding the disposal of accessories subject to wear. To dispose of the appliance, refer to Directive 2012/19/EU.
28. Do not use the appliance if, after a fall, it shows signs of damage on any of its parts. If in doubt, do not use the device, and contact the your country's Mediblink representative company Service Centre.
29. In the case of failure or malfunction of the appliance, switch it off by removing the plug from the mains and contact the Mediblink representative Service Centre. Do not attempt to open, repair or tamper with the appliance.

2-stage nebuliser cup

The effectiveness of aerosol therapy depends on the quality of medicine diffusion through the respiratory tract. The nebuliser cup plays an essential role in creating particles with a diameter that is optimal for therapy, ensuring fast and uniform treatment.

Mediblink 2-stage nebulizer cup produces particles as small as $2,74 \mu\text{m}$ (MMAD), so approximately 86 % of particles achieves lower airways, all the way to lungs alveoli (PORR).

Mediblink 2-stage nebulizer cup also enables the user to set the nebulization rate to two different levels: MIN (0,2 ml/min) and MAX (0,254 ml/min).

The nebulization rate (speed) can be changed by pressing the **BUTTON** (options MIN/MAX) on the top of the nebulization cup (Picture 3).


Table 1: Nebulizer cup nebulization rate settings: MIN / MAX

BUTTON on the top of the cup position	MIN	MAX
NEBULIZATION RATE	$\geq 0.2\text{ml/min}$	$\geq 0.254\text{ml/min}$
NEBULIZATION TIME (5ml of saline solution)	25 min	20 min

The button on the top of the nebulizer cup can be switched when the nebulizer is in use.

2-stage nebulizer cup

Picture 3: Nebulizer cup parts

2 Masks: for adults and children

Soft masks which adapts to your face, minimising drug waste. The masks were designed to ensure correct and efficient therapy for children and adults alike.

Adults should breath through the mouth, as this is the most efficient way to deliver the flow of solution directly into the lower airways, while the children mask should cover both nose and mouth – this is due to the fact that children are not yet able to coordinate their breathing.

Mouthpiece

The mouthpiece can be used alone, without the mask. The mouthpiece conveys the aerosol jet directly into respiratory tract. This minimises drug waste while maximising efficacy.

Preparing and using the appliance

1. Take the appliance out of the package;
2. Insert the medicine and/or saline solution in the nebuliser cup, according to the doses recommended by your physician, without exceeding the maximum level;
3. Close the nebuliser cup and take the air tube and attach it to the compressed air outlet on the nebulizer;
4. Connect the air tube to the nebuliser cup;
5. Connect the nebuliser cup to the accessory needed for the therapy: mask for adults and children, or mouthpiece. In case of doubt, ask your physician about which accessory to use.
6. Connect the power cord plug to the mains;
7. Turn on by pressing the on/off switch and start therapy;
8. Start with nebulization;
9. Carry out treatment preferably seated and in a comfortable position;
10. When the aerosol jet becomes intermittent, interrupt the treatment for a few seconds to allow the suspended drops of solution to fall off the nebuliser cup walls. Resume the treatment session, which will be finished when no more nebulised solution is released from the nebuliser cup;
11. Once the treatment session is completed, switch the appliance off, unplug it, remove the accessory used and clean the appliance and accessories according to the instructions in the Cleaning and Maintenance paragraph.

Cleaning and maintenance

Cleaning and maintenance of the compressor nebulizer

Warning!

Do not wet the device or put it into water or other liquids. For cleaning, use a clean, dry cloth only.

Cleaning and maintenance of the nebuliser cup and accessories

After each use, remove and detach all the nebulizer components and separately flush each part (except air tube) of the nebulizer under running water. **Water temperature must not exceed 50 °C.** Allow all nebulizer parts to air-dry thoroughly.

Do not clean the air tube. If the water droplets accidentally enter the air tube, and the natural wind can not dry the air tube, then only air tube must be connected to the nebulizer, and the nebulizer must be switched on to discharge the water droplets from the air tube. When no droplets are left in the air tube, the nebulizer cup and mask can be attached on the air tube, and the user can start using the nebulizer.

Warning!

Cold disinfect all parts that come in contact with the patient using denatured alcohol. Carefully dry all components and the device before reassembling. Store in a cool, dry place, away from light and heat.

Warning!

Never use benzene, thinners or other inflammable chemical substances for cleaning.

For increased hygienic safety, never use the same accessories on more than one patient. Contact local Mediblink representative to order nebulizer spare parts.

Checking and replacing filters

Periodically check the condition of the filter. The filter was inserted to protect the compressor. Correct filter maintenance prolongs the life of the nebulizer. The frequency with which filter should be replaced depends on the conditions in which the unit is used.

The filter is located on the front of the nebulizer.

To replace the filter:

1. Open the filter cover cap;
2. Remove the filter to be replaced from its housing;
3. Carefully fit the new filter;
4. Fit the filter cover cap back on.

Warning!

Never leave the filter cover cap and/or filters unattended: These are small parts posing the risk of suffocation if swallowed by children.

Troubleshooting guide

Problem	Probable Cause	Solution
Too much noise.	No filter in the unit.	Put a filter in the unit.
The Nebulizer does not run.	Main compressor switch is off.	Turn the unit on.
	Nebulizer was not cleaned after its previous use.	Clean the unit.
	Air tube is bent.	Remove any bend or kink in the air tube.
	The filter is clogged.	Replace the filter.
	No medicine is left.	Add the appropriate amount of medicine prescribed by your physician to the Nebulizer cup.
	The cap is not attached properly.	Reattach the cap.
Water drops form in the air tube.	There is too much medicine in nebulizer cup.	Adjust the amount of medicine in the Nebulizer cup to the appropriate level, connect the air tube to the Air Compressor, and turn the unit on.
	The Nebulizer accessories were washed but not dried.	Dry the Nebulizer accessories until the drops disappear.
NOTE: Please contact your dealer or service company for assistance in case you have any other difficulties.		

Product Specification

Feature	Model / REF	JLN-2321AS / M440
1. Appearance		
2. Device Type	Piston Compressor Nebulizer	
3. Operating Voltage	230 V AC, 50Hz	
4. Power	$\leq 65\text{W}$	
5. Flow	$\geq 8\text{L/min}$	
6. Noise	$\leq 70\text{dB}$	
7. Particle Size	$0.5\text{-}6\mu\text{M}$	
8. MMAD	$< 5\mu\text{M}$	
9. Medication Capability	2~8ml	
10. Working pressure-nebulizer	50kPa-150kPa	
11. Max Pressure	$\geq 230\text{Kpa}$	
12. Working temperature	$10^{\circ}\text{C}\text{-}40^{\circ}\text{C}$	
13. Working humidity	15%-90%RH	
14. Storage temperature	$-25^{\circ}\text{C}\text{-}55^{\circ}\text{C}$	
15. Storage humidity	10%-95%RH (non-condensing)	
16. Atmospheric pressure	700 hPa -1060hPa (working); 500 hPa -1060hPa (storage and transportation)	
17. Weight/Dimension	About 1.4kg / 187mm*142mm*93mm	
18. Life time	5 years	
19. Packaging requirements	Meet the vibration & drop requirements	
20. Disposal requirement	Comply with WEEE directive, it must be disposed of in accordance with the locally applicable regulations, not with domestic waste.	
21. Safety requirements	Comply with EN 60601-1 standard	
22. EMC requirements	Comply with EN 60601-1-2 standard	
23. Performance requirements	Comply with EN ISO 27427 standard	
24. Software	No programmable electronic software is contained	

Product Specification

Information on the operation of the device is provided by the manufacturer and may not be relevant to suspension or high viscosity medicinal products according to EN 13544-1.

EN 60601-1 Medical electrical device — Part 1: General safety requirements

EN 60601-1-2 Medical electrical device — Part 1: General safety requirements

– Safety Standard: Electromagnetic Compatibility – Requirements and Tests

EN 13544-1 Nebulization device — Part 1: Spraying systems and their components.

EMC Information

This section of the instructions contains information relating to the conformity of the product with EN 60601-1-2. The Mediblink Compact M440 nebulizer is a medical device that requires special attention with regard to electromagnetic compatibility and should therefore be installed and used according to the instructions and information provided in the package leaflet. Mobile and portable RF devices (mobile phones, etc.) can affect the operation of the device.

Guidance and manufacturer's declaration – electromagnetic emissions		
The Mediblink Compact Nebulizer M440 (JLN-2321AS) is intended for use in the electromagnetic environment specified below. The customer or the user of the M440 should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The M440 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The M440 is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

EMC Information

Guidance and manufacturer's declaration – electromagnetic immunity			
The M440 is intended for use in the electromagnetic environment specified below. The customer or the user of the M440 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6kV Contact ±8kV Air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %. If ESD interfere with the operation of equipment, counter measurements such as wrist strap, grounding shall be considered.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for Power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1kV differential mode ±2kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 sec	<5% UT for 0.5 cycle 40% UT for 5 cycles 70% UT for 25 cycles <5% UT for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the M440 requires continued operation during power mains interruptions, it is recommended that the M440 be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

EMC Information

Guidance and manufacturer's declaration – electromagnetic immunity – for EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

Guidance and manufacturer's declaration – electromagnetic immunity			
The M440 is intended for use in the electromagnetic environment specified below. The customer or the user of the M440 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 V	Portable and mobile RF communications equipment should be used no closer to any part of the M440, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P} \text{ 80 MHz to 800 MHz}$ $d = 2.3 \sqrt{P} \text{ 800 MHz to 2.5 GHz}$ where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 

EMC Information

Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM – for EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

Recommended separation distances between Portable and mobile RF communications equipment and the M440			
The M440 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the M440 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the M440 as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	50 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1,2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Legend of symbols

	Warning!
	Consult instructions for use Electronic instructions for use: http://www.mediblink.com/f/m440.pdf
	Unique Device Identifier
	Product reference number
	Model Number
	Medical device
	Batch number*

	Device with Type BF applied part
	Class II appliance
	Alternate current
	Keep dry
	Number of products in one packaging
	Manufacturer
	Authorized Representative in EU

*The production date can be read from LOT number [YYYYMM]; first four digits represent the Year and last two digits Month of production. Example: LOT 202503 = March 2025

Conditions for storage and transportation

	Humidity limitation: 10% / 95%
	Keep dry (Store in a cool, dry place)

	Temperature limit: -20 °C / +70 °C
---	------------------------------------

Environmental conditions during use

	Humidity limitation: 10% / 95%
---	--------------------------------

	Temperature limit: +10 °C / +40 °C
---	------------------------------------

This product complies with Directive 2012/19/EU.

The crossed bin symbol on the appliance indicates that the product, at the end of its life, must be disposed of separately from domestic waste, either by taking it to a separate waste disposal site for electric and electronic appliances or by returning it to your dealer when you buy another similar appliance. The user is responsible for taking the appliance to a special waste disposal site at the end of its life. If the disused appliance is collected correctly as separate waste, it can be recycled, treated and disposed of ecologically; this avoids a negative impact on both the environment and health, and contributes towards the recycling of the product's materials. For further information regarding the waste disposal services available, contact your local waste disposal agency or the shop where you bought the appliance.

Warranty

Product: Mediblink Compressor Nebulizer Compact M440

Manufactured for (Importer for EU & Distributor): Mediblink d.o.o.,
Gubčeva cesta 19, 8210 Trebnje, Slovenia; info@mediblink.com; www.mediblink.com

Sellers name, address, signature and stamp*:

Date of extradition / sales*:

*If the invoice is accompanied by this warranty, and if all above information can be seen from the invoice, it is not necessary to fill in this field.

WARRANTY TERMS

Dear customers!

The warranty period is **10 years** and starts from the day of product purchase. The warranty period on spare parts: masks, mouthpiece, tube and filters is **1 year**. In case of product claim, you have to show the invoice. We kindly ask you to save the invoice!

Unfortunately, wrong handling with the device is a reason for 95% of customer complains. You can easily avoid any problem, by getting useful information provided by our special service department. To reach our service department, you can call or send e-mail to Mediblink local distributor.

Before sending the product back to retailer, we kindly ask you to call our service department, to get help about how to use the device to save you with unneeded trips.

The manufacturer guarantees free elimination of all imperfections due to defects in material or manufacturing procedure by repairing or replacing the product. In case that the product can not be repaired or replaced, the customer will get the money refund.

The guarantee is not valid in case of the force majeure, accidents or unexpected events (such as lightning, water, fire etc.), incorrect use or incorrect transport, non-compliance with safety and maintaining regulations or in case of unprofessional product intervention.

Traces of every day product usage (scratches, abrasions) and not subject to claim. The warranty does not eliminate the customer rights, which originate from seller responsibility for product flaws. By accepting the claimed product by the service department, the service department does not take responsibility for loss of saved data or settings on the product. All product repairs, which are performed out of product warranty period, have to be paid by customer by prior notice.

The manufacturer guarantees the product quality and flawless product operation in the warranty period, which starts with the day of product purchase. If the product can not be repaired in 45 days, the product will be replaced with a new one. In case that the product can not be replaced, the money will be refunded to the customer.

In case of product claim, call or send e-mail to Mediblink local distributor.

Warranty

Replacement parts can be ordered at www.mediblink.com.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Manufactured for
(Importer for EU & Distributor):
Mediblink d.o.o.
Gubčeva cesta 19
8210 Trebnje
Slovenia
info@mediblink.com
www.mediblink.com

Shenzhen Homed Medical Device CO., Ltd.
Room 1601-1614, A Building, Jinxiu III, No. 85 Hudi Pai,
Songxuan Community, Guanhu Street, Longhua District,
Shenzhen 518109, P. R. China
E-mail: info@homedgroup.com

Shanghai International Holding Corp. GmbH (Europe)
Eiffestrasse 80, 20537 Hamburg, Germany
E-mail: shholding@hotmail.com

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